



officine
orobiche

quantity 14

off No.

P18/2683 pos 2 of 2 rev 0

model CVM20.JL.MO.M4.WP.XX

tag

s1
s2
s3
s4
s5

s6
s7
s8
s9
s10

s11
s12
s13
s14
s15

s16
s17
s18
s19
s20

01
02
03
04
05

body mat. AISI 316L
connections 3/4" ANSI 150 RF 125-250 AARH
float MONEL
attraction sleev Nickel
flow indicator YES

21
22
23
24
25

body rating
max-press. 10 barg
max-temp. 100 °C
fluid H2O

06
07
08

housing ALUMINIUM WP IP 66
AREA .

26
27
28

nominal flow 2 m3/h
range 0.4---1 m3/h
set (min) 0,6 m3/h

09
10
11
12

electric connection 1/2 " NPTF
switch groups 1
microswitch M4
contact SPDT (125 V DC 0,5 A)

29
30
31
32

diff. 0,6 m3/h
normal flow 2 m3/h
adjust.(min) on site
scale CM 29 (mc/h)

13
14
15
16
17

connector
cable gland
fins
label
paint
ALUMINIUM

33
34
35
36
37

nace
NO

18
19
20


tube
internal parts
antisteam
MONEL
MONEL

38
39
40

mat_cone
MONEL

v1
v2
v3
v4
v5

n1
n2
n3
n4



gen-18

issued UV.....date

approved UT.....date