



**officine
orobiche**

quantity 1 off No. P18/2930 pos 2 of 2 rev 0

model PLD100.J.2.M2.WP

tag

s1
s2
s3
s4
s5

s6
s7
s8
s9
s10

s11
s12
s13
s14
s15

s16
s17
s18
s19
s20

01
02
03
04
05

body mat.
connections
paddle
attraction sleeve
flow indicator

AISI 316
4" ansi 150 RF 125/250 aarh
AISI 316
AISI 430

21
22
23
24
25

fluid
max-press.
max-temp.
hydraulic test

H2O
10 barg
100 °C
15 barg

06
07
08

housing ALUMINIUM WP IP 66
.
AREA .

26
27
28

nominal flow
normal flow
set (min)
diff.

60 mc/h
60 mc/h
12 - 30 mc/h
30 mc/h
23 mc/h on site

09
10
11
12

electric connection
switch groups
microswitch
contact

1/2 " NPTF
1
M2
SPDT (220 V AC 10 A)

29
30
31
32

adjust.(min)

13
14
15
16
17

connector
cable gland
fins
label
paint

.
.
.
ALUMINIUM
.

33
34
35
36
37

nace

no

18
19
20


tube
internal parts
antisteam

AISI 316
AISI 316
.

38
39
40

v1
v2
v3
v4
v5

n1
n2
n3
n4



gen-2018

issued UV.....date