




officine
orobiche

quantity 2

off No.

P18/3329 pos 1 of 1 rev 0

model 72D.A.2.M2.EP.XX

tag		fluid	Density Kg/dm3	temp. es / prog °C	press. es / prog barg	switch point H - L	
s1	425LSLL026	LIQUID	1	20 / 40	1 / 1	. - - -	
s2	425LSH025	LIQUID	1	20 / 40	1 / 1	. - - -	
s3							
s4							
s5							
01	body mat.	CARBON STEEL			21	body rating	####
02	connections pos.	SIDE - SIDE 356 mm			22	hydraulic test	####
03	connections type	1"ANSI 150 RF 125-250 AARH			23		
04	drain	1 " NPTF			24	float	AISI 316
05	vent	.			25		
06	housing	ALUMINIUM ATEX II 2 G Ex d IIC T6 o T5 Gb			26		
07		II 1/2 D Ex ta/tb IIIC T 85°C o T100°C Da/Db			27		
08	AREA	Body ATEX II 2 GDc Tx			28		
09	electric connection	1/2 " NPTF			29		
10	switch groups	N. 2			30		
11	microswitch	M2			31		
12	contact	SPDT (24 V DC 1 A)			32	nace	NO
	connector	.					
13	cable gland	.			33	specific customer	.
14	fins	.			34	specific customer	.
15	paint	stc 088 (silic.)			35		
16	bolts	A193.B7 GALVANIZED			36	label	ALUMINIUM
17	attraction sleev	AISI 430			37		
18	tube	AISI 316			38		
19	internal parts	AISI 316			39		
20	antisteam	NO			40		
v1	XX: Center to center distance= 400 mm						
v2							
v3							
v4							
v5							
n1							
n2							
n3							
n4							
n5							
A	.						

issued UV.....date

lug-18

issued UV.....date

lug-18